




TapeTx Medical Taping Basic Skills

1. Concept



2. Purpose

- 1) Provide the strength to close the mouth
- 2) Provide the fixation of the oral immobilizer to prevent the backward movement of the jaw
- 3) No glue area to prevent lip injury
- 4) No glue area to provide using safety
- 5) No glue end point to provide easy and painless tape removal

3. Using Timing

- 1) Attach before sleeping
- 2) Detach after waking up
- 3) Reusing is not suggested due to sanitary consideration

4. Attaching Procedure

- 1) Make a cut around 5-7 cm (meet the width of lips)
- 2) Stretch the width of no glue area till enough to protect the lips
- 3) Rotate the tape to cover the no glue area on the lips
- 4) Attach the 1st glue area above the upper lip
- 5) Attach the 2nd glue area below the lower lip

5. Detaching Procedure

- 1) Insert the finger into one end opening of the no glue area
- 2) Pull the end opening to the other side of the tape to detach the tape

6. Precautions

- 1) TapeTx Medical Taping should be applied after the diagnosis and evaluation of the medical doctor.
- 2) The patented no glue design and the glue of the TapeTx Medical Tape is different from the ordinary tapes. Do not use ordinary tapes to do the taping.
- 3) Do keep the no glue area parallel to the lips to protect them.
- 4) The oral immobilizer might cause biting discomfort at the beginning. Please let your medical doctor know and evaluate when you feel it.
- 5) People of severe nasal and inspiration track obstruction are not suggested to do



TapeTx Medical Taping.

- 6) People of stroke or upper limb mobilization disability are not suggested to do TapeTx Medical Taping.
- 7) Children less than 4 years old are not suggested to do TapeTx Medical Taping.
- 8) People with skin irritation or breathing difficulty after taping are not suggested to do TapeTx Medical Taping.

OMNIBUD